

QUICK CLINICAL WORKFLOW GUIDE

A simple, infection-based roadmap to guide bacterial testing and treatment decisions

NEW PERIO PATIENT (Diagnosed Today)

Who this applies to:

New patients diagnosed with periodontal disease based on:

- Pocketing
- Bleeding
- Inflammation
- Radiographic bone loss

What to Do

- At diagnosis, collect a paper-point sample from the deepest bleeding pockets
- Send the sample to MicrobeLink Dx®
- Results available in ~72 hours
- If the patient is ready to begin care, SRP may begin immediately
- Scaling and root planning is always performed for perio patients

When to Consider Systemic Antibiotics

Antibiotics are optional, not automatic. Consider if:

- Infection has been present long-term
- Bacterial load is high (especially red complex)
- Patient has high systemic risk (heart disease, diabetes, autoimmune, etc.)

 If infection is not severe:

- Complete SRP without antibiotics
- Re-test in 3–6 months
- Escalate treatment later if needed

Follow-Up

- Re-test in 3–6 months, depending on severity

Verbal Skills

“You have an active infection. This bacterial test tells us how severe the infection is and how aggressively we need to treat it”

REFRACTORY PERIO PATIENT (Not Responding to SRP)

Who this applies to:

Patients who have completed SRP without bacterial testing and still present with:

- Bleeding
- Inflammation
- Persistent pocketing

What to Do

- Collect a paper-point sample from the deepest bleeding pockets
- Send the sample to MicrobeLink Dx®
- Results available in ~72 hours

When to Consider Systemic Antibiotics

Refractory cases are more likely to require antibiotics. Consider if:

- Mechanical therapy alone has failed
- Report shows high bacterial load or aggressive pathogens
- Patient has systemic risk factors

 Antibiotics are most effective when combined with biofilm disruption

Follow-Up

- Re-test 3–6 months

Verbal Skills

"We've tried scaling and root planning, but the infection is still there. This test shows us why and helps us be more aggressive where needed."

6-MONTH RECALL PATIENT (NOT A PROPHY – "CONVERSION CASE")

Who this applies to:

Patients presenting for prophy recall with:

- Bleeding
- 5 mm pockets
- Inflammation
- New or worsening medical history

What to Do

- Recognize this is not a prophy
- Collect a paper-point sample from the deepest bleeding pockets
- Send the sample to MicrobeLink Dx®

Treatment Approach

- Treat like a new perio diagnosis
- SRP is indicated
- Adjunctive therapy based on bacterial results

Antibiotics

- Consider only for severe infection or systemic risk

Follow-Up

- Re-test in 3–6 months

Verbal Skills

"This isn't just inflammation. This test shows us what bacteria are keeping your gums from staying healthy."

IMPLANT CANDIDATE (Before Placement)

Who this applies to:

Patients scheduled for implant placement, especially with:

- History of perio
- Bleeding or inflammation
- Bone loss
- Systemic health concerns

What to Do

- Before implant placement, collect a paper-point sample from the implant quadrant
- Send the sample to MicrobeLink Dx®

Treatment Approach

- Treat infection the same way as a perio patient
- If SRP is not indicated, disrupt biofilm prior to antibiotics
- Goal: reduce bacterial load before surgery

Antibiotics

- Consider if bacterial load is moderate–high
- Combine with biofilm disruption (not antibiotics alone)

Follow-Up

- Re-test before implant placement to confirm infection control

Verbal Skills

"We never want to place an implant into infection. This test helps give you the best possible outcome."

FAILING IMPLANT / PERI-IMPLANTITIS

Who this applies to:

Implants with:

- Bleeding
- Mucositis

- Inflammation
- Bone changes

What to Do

- At the first sign of tissue changes, collect a site-specific paper-point sample
- Send the sample to MicrobeLink Dx®

Treatment Approach

- Treat similarly to refractory perio cases
- Disrupt biofilm around the implant
- Use results to guide aggressiveness of care

Antibiotics

Often more appropriate in failing implant cases if:

- Bacterial load is high
- Aggressive pathogens are present

Follow-Up

- Re-test in ~3 months

Verbal Skills

"Implants can get infections just like teeth. This bacterial test tells us how severe the infection is and how aggressively we need to treat it"

SYSTEMICALLY HIGH-RISK PATIENT

Who this applies to:

Patients with conditions such as:

- Heart disease
- Diabetes
- Pregnancy or pregnancy planning
- Alzheimer's risk
- Autoimmune or immunocompromised conditions
- Joint replacement

What to Do

- At the first sign of bleeding or inflammation, collect a paper-point sample
- Send the sample to MicrobeLink Dx®

Treatment Approach

- Treat proactively, even with mild pocketing
- Adjust treatment intensity based on bacterial severity and medical risk

Antibiotics

Consider if:

- Bacterial load is high
- Systemic condition makes infection unsafe

Follow-Up

- Re-test every 3–6 months

Verbal Skills

“Because of your medical history, controlling oral infection is especially important. This test helps us protect your overall health.”

MicrobeLink Dx® | Paper-point DNA-PCR testing | The Science Behind the Smile
www.microbelinkdx.com